HARTMANN

Tensoval®

mobil

Instructions for use
Instrucciones de uso
Manual de utilização
Οδηγίες χρήσεως
Kullanım Kılavuzu
تعليمات للاستعمال

 Please read this introduction carefully before use as correct measurement can only be carried out following these instructions. Should you have any further questions on our product or on measuring blood pressure, please contact our Customer Service (a list of addresses is given in the warranty document). On the Internet, you can visit our homepage on www.hartmann.info.

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Congratulations on your purchase of this HARTMANN quality product. Tensoval mobil is a fully automatic blood pressure monitor which uses intelligent measuring technology (Fuzzy Logic) to take rapid and reliable measurement of systolic and diastolic blood pressure as well as pulse rate using an oscillometric measuring method.

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
1. General information on blood pressure

The human heart beats approx. 60 to 80 times per minute. This causes blood to be pumped into the arterial vascular system and supplies the body with oxygen and necessary nutrients. A constant pressure known as blood pressure is required to enable the blood to flow into the smallest of blood vessels. In measuring blood pressure, two values are determined. The **systolic value** is the maximum pressure as the heart beats. **Diastolic value** is the

minimum pressure between two heart-beats. Blood pressure is expressed as, e.g. 120/80 in millimetres of a mercury column (mmHg).

Blood pressure changes constantly in each individual and thus creates the requirements necessary for the functional capacity of the body. Thus fluctuations in blood pressure are perfectly normal. However, if blood pressure values are constantly elevated at rest, **hypertension** or high blood pressure is present.

The World Health Organisation (WHO) has set out the following limit values:

Assessment	Systolic pressure	Diastolic pressure
optimal	up to 120 mmHg	up to 80 mmHg
normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg

High blood pressure is one of the commonest causes of disability and death. And yet many people do not know they are suffering from high blood pressure. This is because the problem is only noticed when the conditions it leads to are at an advanced stage. These may be life-threatening, such as heart attack, stroke or kidney damage. Only regular blood pressure monitoring enables high blood pressure to be detected early. You have decided in favour of **blood pressure self-measurement** and this will contribute to your health care. This is the only way to ensure efficient protection against the dangerous results which may affect your heart and circulation.

2. Important notes for self measurement

- Self measurement of blood pressure does not replace either treatment or **regular check-ups by your doctor**. The drug dose prescribed by your doctor may thus on no account be altered without consultation.
- Self-measurement of blood pressure must be conducted regularly and over a long period of time. With the assistance of the **blood pressure values recorded regularly by you in your blood pressure diary**, high blood pressure can be detected early and your doctor can initiate optimal treatment. Individual readings are situation-related and thus are not useful.
- In order to compare values, it is important that you **always** take your blood pressure **at rest**. So you should relax for 5 minutes before taking a measurement.
- Even slight changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to **fluctuations in blood pressure**. This is why your doctor and pharmacist often obtain different readings.
- Measure your blood pressure regularly in the morning and evening **at about the same time**, as blood pressure fluctuates during the course of the day.
- You should **not move or talk** during the measurement process.
- Always measure blood pressure on the same arm, and hold your wrist still throughout the entire measuring process **at the level of your heart**. If the position of the wrist is too high or too low, the measurements could be incorrect.
- Tensoval mobilE is intended for a **wrist circumference between 12.5 and 20.5 cm**. Outside these limits correct readings may no longer be guaranteed.
- When taking repeated measurements, **relax for at least one minute** between measurements as otherwise false results are obtained.
- In severe **cardiac arrhythmias** (cases

of disturbed heart rhythm), measurements should only be taken in consultation with a doctor. Due to the oscillometric measuring method, measurements taken with this device may in some instances be false or fail to produce a reading (Err).

- In circulatory disorders or severe **arteriosclerosis**, it is not possible to take a reliable reading of blood pressure on the wrist. In such cases, the upper arm should be used for measurement.

- The monitoring of blood pressure is extremely important in **pregnancy**, as the pressure may be changed by the pregnancy. The results should, however, only be interpreted in consultation with the doctor.

3. Control displays and symbols

- ♥ flashes when the device is measuring and the pulse is being taken.
- 🔋 change batteries
- Err measuring error, cf. Chap. 12
- ⬆ appears during inflation
- ⬇ appears during automatic checking

M1 displays the saved reading for person 1

M2 displays the saved reading for person 2


protection from electric shock (type BF)

pay attention to operating instructions

4. Current supply

- Open the battery cover on the left side of the device (see fig.) Insert batteries (see chap. 13, Technical data), ensuring correct polarity ("+" and "-"). Close the battery lid.

- We recommend the use of high-quality batteries as other batteries or accumulators could affect the precision

of measurement. Never mix old and new batteries or batteries made by different manufacturers.

- If you do not intend to use the device for a longer period of time, please remove the batteries.
- In the interests of environmental protection, used batteries should not be disposed of in household waste. Please observe valid regulations when disposing of the batteries or use municipal collecting bins.


5. Setting date and time

- After inserting the batteries the time function automatically appears. The flashing year key is then displayed. The year can be changed using the M1 (+) and M2 (-) buttons. Save the **year** by pressing the START/STOP button.
- Next the **month** is saved. The number on the right will flash. Use the buttons as above for setting the year. Repeat procedure for saving the **day, hour and minute**. Each time you change the batteries, the date and time must be reset.

6. Applying the blood pressure watch

- Measurement should be taken on the arm which usually has the higher blood pressure value. Remove your watch or other wrist jewellery.
- Place the Tensoval mobil on the

inside of the wrist. The edge of the cuff should lie about 1 cm from the wrist fold. Close the cuff and pull as firmly as possible (see fig.) and close the Velcro fastener securely.


- The cuff is designed for **wrist circumferences of 12.5 to 20.5 cm**. Check your wrist circumference with the assistance of the red strip at the edge of the cuff. The point should be within the marking strip.
- The Velcro fastener tip may be used to tuck away the extra length of cuff.

7. Measuring blood pressure


- We recommend that blood pressure should be taken in the sitting position. Ensure that the **wrist is at the level of the heart, relaxed on the chest**, and is not moved during measurement. The other hand may be used to give the

elbow support in order to relax the muscles (see fig.).

and holding it down until the desired cuff pressure is reached. This should lie about 30 mmHg above the systolic value.

- **Important: do not move or talk throughout the entire measuring procedure.**
- As pressure in the cuff decreases, the heart symbol and the falling cuff pressure are displayed.
- A **beep** indicates the end of measurement. The **systolic and diastolic readings** appear simultaneously on the display as well as the pulse rate beneath (see fig.).


The **time** appears above the reading and M1 or M2 is displayed on the left. **M1** represents the reading for the **first person** and **M2** for the **second**. As long as the reading is displayed, you can apply the values to the respective person by pressing M1 or M2. Otherwise, if they are not applied the reading is automatically saved for the person on the display.

- Do not switch the device on until it has been attached to the wrist, as the cuff may otherwise be damaged by excess pressure.
- Press the **START/STOP button**. The appearance of all display segments followed by a flashing downward-facing arrow, shows that the device is checking itself automatically and is ready for use. Afterwards, **automatic inflation** begins up to about 190 mmHg. If this inflated pressure is insufficient or measurement is interrupted, the device will pump in 30 mmHg stages up to the appropriate higher pressure.
- If as an individual you require a higher inflation pressure, you can circumvent the subsequent inflation, by pressing the START/STOP button again just after the beginning of the inflation process

- To switch the device off simply press the START/STOP button, otherwise the device will switch itself off automatically
- If you wish to stop measurement for any reason, simply press the START/STOP button. The inflation or measuring process is interrupted and an automatic fall in pressure occurs.

8. Memory function

- To call up stored readings, press the **Memory button** with the device switched off. Press M1 for the first person's recorded values and M2 for the second person. The corresponding symbol M1 or M2 will be displayed. First, the average value of all saved data will be shown for this person. An "A" appears on the display and the number of readings, from which the average value was calculated is displayed in the top right-hand corner (see fig.).

Press the memory button to access memory position No. 1

- Tensoval mobil stores up to **30 measurements**. The latest reading is always in position No. 1, all older readings are pushed back one position. If all memory positions are occupied, the oldest is deleted.
- To call up all stored readings one after another, **press the Memory button repeatedly**.
- The measurement and the relevant number of the memory position appear when you access stored values. The memory position, date and time appear in 2 – 3 second intervals.
- You can interrupt storage function at any time by pressing the START/STOP button. Otherwise the device switches itself off automatically after a few seconds.
- Should the power supply fail, e.g. when changing the batteries, the stored values will still be available.

Deleting stored values

You can delete all stored data in M1 or in M2 for the respective person. To do this, press the corresponding memory button. The average value will then appear on the display. Press down on the memory button and after four seconds the display will flash. After eight seconds all data for the respective person will be deleted and only M1 or M2 will be displayed. To avoid deleting any data, release the button before the end of the eight-second time limit.

9. Maintenance

- Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- This device consists of high quality precision components. Do not drop the device and protect it against hardknocks.
- The moulded cuff should not be folded or overstretched.
- Never open the device. Repairs may only be carried out by qualified individuals.
- Only use a soft, moistened cloth to clean the device. Do not use detergents or solvents.
- The cuff can be cleaned carefully with a moist cloth and mild, pH-neutral soap. Do not completely immerse the cuff in water.

10. Checking of measurement function

We recommend that measurement function be checked at intervals of two years in the case of professional users, such as pharmacies, medical practices and clinics. You should also observe the national regulations, such as, in Germany, the Decree on Sales of Medical Products (Medizinprodukte-Betreiberverordnung). Measurement function checks can be carried out either by Paul HARTMANN AG (see Service address), by competent authorities, or authorised maintenance providers against compensation.

Instructions for Calibration Mode:




Remove the batteries in order to switch to calibration mode. Hold down the START/STOP button and then insert the batteries. Shortly after the button is released, two zeros, one above the other, appear on the display.

Instructions on checking measurement function will be supplied on request to competent authorities, or authorised maintenance providers.

11. Warranty

This product comes with a **3-year warranty**. Details of warranty conditions and contact addresses can be found in the separate warranty document.

12. Error displays

Error which has occurred	Possible causes
Device will not switch on	No batteries, they have been inserted incorrectly or are dying.
	Moving or speaking during the measurement.
	Cuff was incorrectly or loosely applied.
	Please repeat reading.
	Pressure in cuff exceeds 330 mmHg. Therefore, pressure will be automatically released.
	Batteries are almost dead. Only a few more measurements are possible.

Error which has occurred	Possible causes
	Batteries are dying. Insert new batteries.
Implausible reading	The blood pressure watch has been positioned wrongly. Wrist circumference not within permitted range. Wrist not held at heart level. Moving, talking or excitement during measurement. Deep breathing during measurement process. Lack of relaxation before measurement. Stimulants taken just before measurement.

Switch the device off if an error symbol appears. Check the possible causes and note the instructions on self-measurement in Chapter 2. Relax for a minute and take the measurement again. Do not move or talk during the measurement.

13. Technical specifications

Measuring method:	oscillometric
Display range:	0 – 300 mmHg
Measuring range:	Systole (SYS): 50 – 250 mmHg Diastole (DIA): 40 – 180 mmHg Pulse: 40 – 160 beats/ minute
Technical accuracy of measurement:	Cuff pressure: ± 3 mmHg Pulse: ± 5 % of displayed pulse frequency
Clinical accuracy of measurement:	complies with the EN 1060 standard, part 3.
Power supply:	2 x 1.5 V alkaline Mangan Micro (AAA/LR03) batteries
Battery capacity:	> 2,000 measurements
Inflation pressure:	approx. 190 mmHg
Automatic cut-out:	3 minutes after end of measuring process.
Cuff circumference:	12.5 – 20.5 cm
Pressure release valve:	electronically regulated and pulse-controlled linear valve
Storage capacity:	2 x 30 readings and mean value
Operating conditions:	Ambient temperature: $+10$ °C to $+40$ °C Relative humidity: 15 – 90 %
Storage/transport conditions:	Ambient temperature: -20 °C to $+50$ °C Relative humidity: 15 – 90 %
Serial number:	see battery compartment

14. Legal requirements and guidelines

Tensoval mobil complies with the requirements of the EC directive, 93/42/EEC on Medical products and bears CE-marking.

The device complies with the EN 1060 standard relating to non-invasive blood pressure monitors Part 1/1995: General requirements and Part 3/1997: Additional requirements for electromechanical blood pressure measuring systems.

Manufacturer:
PAUL HARTMANN AG, 89522 Heidenheim
Germany

15. Spare parts

To ensure measurement accuracy, only use original HARTMANN spare parts which may be obtained from your pharmacist or specialist medical supplier.

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DE – PAUL HARTMANN AG · 89522 Heidenheim
EG – PAUL HARTMANN Egypt S.A.E. · Giza, 12515
ES – Laboratorios UNITEX-HARTMANN S.A. · 08302 Mataró
GB – PAUL HARTMANN Ltd. · Heywood/Lancashire OL10 2TT
HK – PAUL HARTMANN Asia-Pacific Ltd. · Hong Kong
SG – PAUL HARTMANN Pte. Ltd. · Singapore 609 143